






**amtelco**




R&D Software Department  
 4800 Curtin Drive, McFarland, WI USA 53558  
 www.amtelco.com






















## New Features in Amtelco Secure Messages Contact Web 7.0

All rights reserved © May 2025

FEATURE	DESCRIPTION
<p>Accounts Links in Navigation Menu            7.0.9144.1</p>	<p>The Amstelco Secure Messages navigation menu displays the accounts to which you are registered.</p>  <p>Links for the Amstelco Secure Messages accounts that you are registered to are displayed at the bottom of the navigation menu.</p> <p>To switch to a different Amstelco Secure Messages account, click the name of the account.</p>
<p>Role-Based Messaging            7.0.8892.2  <i>optional feature</i></p>	<p>The optional Role-Based Messaging feature has been implemented.</p> <p>Role-Based Messaging synchronizes Intelligent Series (IS) OnCall Scheduling with Amstelco Secure Messages. Role-Based Messaging to enable Amstelco Secure Messages users to retrieve assigned resources for IS on-call schedules, on-call shifts, and on-call roles and to send messages to everyone filling a role.</p> <p><b>Note:</b> This feature requires Amstelco Secure Messages Admin Web 7.0 or later and Intelligent Series 5.7 or later.</p>

FEATURE	DESCRIPTION
	
	<p data-bbox="548 898 1437 976">The Contacts Screen now displays roles as well as contacts. Roles added to a Contact Circle receive messages sent to the Contact Circle.</p> 
	<p data-bbox="548 1522 1372 1591">The Roles Screen is used to view the roles that you are currently filling and to view and send messages to contacts filling roles.</p>

FEATURE	DESCRIPTION
	 <p>The Compose Message window contains the Send As menu, allowing users to choose to send a message in their own name or to send a message as a role they are filling. Messages sent to or from a role stay with the role, not the user, allowing message threads to carry over from one shift to the next for seamless coverage.</p>
<p>Away Message duration error prompt 7.0.8892.2</p>	<p>When setting an Away Message in the Settings page, an error prompt is displayed if the amount of time entered exceeds the duration limit.</p>  <p>The error prompt displays the maximum allowed Away Message time configured by your administrator.</p>
<p>Amtelco Secure Messages rebranding 7.0.8641.4</p>	<p>As part of the Amtelco rebrand, the miSecureMessages Contact Web has been renamed to the Amtelco Secure Messages Contact Web and features a new logo and theme color.</p> 
<p>Thread View pagination 7.0.8641.4</p>	<p>If there are more than 25 message threads in the Thread View on the Inbox page, the <u>More Threads</u> hyperlink is displayed at the bottom of the Thread View.</p>

FEATURE	DESCRIPTION						
	<div data-bbox="548 268 1450 655"><table border="1"><tbody><tr><td data-bbox="565 283 1201 346"><input type="checkbox"/>   Erica Sorenson Upcoming surgery</td><td data-bbox="1209 283 1442 310">03/25/2024 4:51:39 PM</td></tr><tr><td data-bbox="565 367 1201 430"><input type="checkbox"/>   Amy Thompson, Alan Caspar Patient update</td><td data-bbox="1209 367 1442 394">03/25/2024 4:43:28 PM</td></tr><tr><td data-bbox="565 451 1201 514"><input type="checkbox"/>   Alan Caspar Lab results</td><td data-bbox="1209 451 1442 514">03/14/2024 11:19:41 AM </td></tr></tbody></table><div data-bbox="906 562 1084 615" style="text-align: center; background-color: #e0f2f7; padding: 5px;"><a href="#">More Threads...</a></div></div> <p data-bbox="548 678 1344 751">Click the <u>More Threads</u> hyperlink to display the next group of message threads in the Thread View.</p>	<input type="checkbox"/>   Erica Sorenson Upcoming surgery	03/25/2024 4:51:39 PM	<input type="checkbox"/>   Amy Thompson, Alan Caspar Patient update	03/25/2024 4:43:28 PM	<input type="checkbox"/>   Alan Caspar Lab results	03/14/2024 11:19:41 AM 
<input type="checkbox"/>   Erica Sorenson Upcoming surgery	03/25/2024 4:51:39 PM						
<input type="checkbox"/>   Amy Thompson, Alan Caspar Patient update	03/25/2024 4:43:28 PM						
<input type="checkbox"/>   Alan Caspar Lab results	03/14/2024 11:19:41 AM 						

**Requirements:**

- Amtelco Secure Messages Server 7.0 or later

**Role-Based Messaging Requirements:**

- Amtelco Secure Messages Role-Based Messaging licenses
- Amtelco Secure Messages Admin Web 7.0.8641.3 or later
- Intelligent Series Supervisor 5.7.8788.3 or later